

# Training on the Agilent Open Access GC/MS

---

*Questions?*

Contact Dr. Allis Chien  
allis@stanford.edu  
650-723-0710



STANFORD UNIVERSITY MASS SPECTROMETRY  
333 CAMPUS DR., MUDD 175  
STANFORD, CA 94305-5080

# Table of Contents

---

- ❑ Preparations – sample & instrument (Slides 3-6)
- ❑ Methods (Slides 7-12)
- ❑ Data Acquisition (Slides 13-25)
- ❑ Data Analysis (Slides 26-31)
- ❑ Peak Integration (Slides 32-34)
- ❑ Data Reporting (Slides 35-39)
- ❑ Contact Info (Slide 40)

# Analytes suitable for GC/MS

---

- In general, GC/MS is used to analyze non-polar compounds soluble in non-polar solvents such as the following:
  - Hexane, Ethyl Acetate, Methanol, Acetonitrile, Acetone
  - High-boiling solvents should be used only with appropriate methods
  - \*\*Pyridine and Chlorinated Solvents MUST be avoided.
  
- The compound(s) to be analyzed should be both volatile and chemically stable at the operation temperature of the injector (280° Celsius)
  - Thermally labile compounds may be pyrolyzed at the max temperature of injector and/or column (320° Celsius)
  
- Water solutions, and organic solutions at pH lower than 2 and greater than 8 are NOT allowed for GC/MS analysis, as they may damage the column

# Sample Preparation


---

- Concentration of the analyte: 1 ng/uL
- The solvent used MUST completely dissolve your compound – no cloudy solutions or precipitate!
- Use 2 mL autosampler vials: Agilent PN 5182-0553
  - Recommended minimum sample volume is 500 uL
- For limited volume samples, use 300 uL inserts: MicroSolv PN 9502S-02CP
  - Minimum sample volume using inserts is 30 uL

# Instrument Prep Checklist

---

- ❑ Make sure the solvent vials have enough solvent as indicated by the blue line on the vial.
- ❑ Solvent A: turret position 11  
Acetonitrile, b.p. 82° C
- ❑ Solvent B: turret position 8  
Ethyl Acetate, b.p. 77° C
- ❑ Waste: positions W, A5  
vials present and empty




Solvents are stored under the fume hood; glass pipets are attached to the bottles

# GC-MS Column

---

- All methods use a capillary column with 3  $\mu\text{L}$  of stationary phase HP-5, (5%-phenyl)-methylpolysiloxane which is equivalent to USP Phase G27
- Column dimensions: length 30 m, ID 250  $\mu\text{m}$ , film thickness 25  $\mu\text{m}$
- More application information is available online at

<http://www.chem.agilent.com/en-US/Products/consumables/columns/gcandgc-ms/jwhp-5/pages/gp42676.aspx>



# Methods

Open Access GC/MS



# GC/MS Methods

---

- General survey method
  - **SUMS\_Standard** (1 ng/uL) (methanol, hexane)
- Focus on compounds by boiling point
  - **SUMS\_HighBP** (heptane, iso-octane, toluene)
  - **SUMS\_MedBP** (ethyl acetate, acetonitrile, acetone)
  - **SUMS\_LowBP** (methanol, pentane, hexane)
- Higher sensitivity methods – same temperature program as Standard method
  - **SUMS\_StdMedSens** (0.5 ng/uL)
  - **SUMS\_StdHighSens** (0.1 ng/uL)

- 
- *Method\_Name* (recommended solvents or concentrations)
  - *All methods inject 1 uL of sample*



# Standard

## Methods

---

- SUMS\_Standard                      Split ratio 1:100
- SUMS\_StdMedSens                Split ratio 1:50
- SUMS\_StdHighSens               Split ratio 1:25
  
- Oven Temperature
- Start:                            35° C, hold for 3.75 min
- Ramp:                            20° C/min from 35° C to 320° C
- End:                              320° C, hold for 7 min
  
- Injector Temperature 280° C

# High Boiling Point Methods

---

- SUMS\_HighBP Split ratio 1:100
  
- Oven Temperature
  - Start: 75° C, hold for 2 min
  - 1<sup>st</sup> Ramp: 40° C/min from 75° C to 120° C
  - Isotherm: 120° C, hold for 5 min
  - 2<sup>nd</sup> Ramp: 20° C/min from 120° C to 320° C
  - End: 320° C, hold for 7 min
  
- Injector Temperature 280° C

# Medium Boiling Point Methods

---

- SUMS\_MedBP                                  Split ratio 1:100
  
- Oven Temperature
  - Start:                                  50° C, hold for 2.00 min
  - 1<sup>st</sup> Ramp:                              30° C/min from 75° C to 100° C
  - Isotherm:                              100° C, hold for 4 min
  - 2<sup>nd</sup> Ramp:                              20° C/min from 100 C to 225 C
  - Isotherm:                              225° C, hold for 1.7 min
  - 3<sup>rd</sup> Ramp:                              40° C/min from 225° C to 320° C
  - End:                                      320° C, hold for 7 min
  
- Injector Temperature 230° C

# Low Boiling Point Methods

---

- SUMS\_LowBP Split ratio 1:100
  
- Oven Temperature
  - Start: 35° C, hold for 5 min
  - 1st Ramp: 20° C/min from 35° C to 180° C
  - Isotherm: 180° C, hold for 2.25 min
  - 2nd Ramp: 40° C/min from 180° C to 320° C
  - End: 320° C, hold for 7 min
  
- Injector Temperature 180° C

# Data Acquisition

---

Open Access GC/MS

# Queuing an Individual Sample: Start with "Queue Mode Enabled"

---

Click the "Queue Mode Enabled" button at the bottom of the screen




The "Secured System Control" Window will appear

Choose

1) Run Individual Sample

2) Click "OK"




# Individual Sample: Log in, choose a "method"

---


The "Enter Login Name"  
Window will appear

- 1) Enter your "e-mail" ID
- 2) Enter your "Password"
- 3) Click "OK"



The "Select Analysis Method"  
Window will appear

- 1) Choose one method
- 2) Click "OK"



# Individual Sample: Submit the Sample

---

The "MDS ChemStation " Window will appear

Display the Vial position


**Place your vial in the position #**  
**Do not overwrite the position**

1) Click "OK"



The "Sample Information" Window will appear

- 1) Enter your "text" Sample Title
- 2) Enter your "text" Information
- 3) Click "OK"





# Individual Sample: Submit the Sample




---

The "MSD ChemStation " window will remind you of the vial position:

Your Vial should be in position #  
Be sure your vial is in position #

1) Click "OK"

If you need to submit another sample:  
Click "Yes", re-start from slide #13  
(no login needed); otherwise Click "No"



# Individual Sample: Verify Sample Submission

After you submit the sample request and your sample vial has been placed in the assigned autosampler position:

- Click on "Report Manager" and verify that your sample posted

Report Manager - "CHEMSTATION"

File Options Help

Start Pause Resume Shutdown

Posted 25 qmeth "37","One Sample"," Single run","sums" on P37669.d

# of Requests 1

System/Command	Data File	Method	Originator
RPT1 on \\CHEMSTATION\ using MSDA			
RPT2 on \\CHEMSTATION\ using MSDA2			
MS 1 on \\CHEMSTATION\ using MSTSP1	25 qmeth "37","One Sample"," P37669.d	SUMS_Standard.M	CHEMSTATION

Vial Position

Filename

Method

- Important! Write the Filename in your lab notebook

# Individual Sample: Important reminders

---

- Never close windows on the instrument control computer! Only minimize them
  - Closing/exiting Chemstation windows will stall the queue and aggravate your colleagues
  
- All queued analyses are posted in the "Report Manager"
  - Never remove a vial in a position which is still listed in the "Report Manager"
  - If a vial position is not listed in the "Report Manager", a vial present in that position may be replaced with your sample vial

# Queuing a Batch of Samples: Start with "Queue Mode Enabled"

---


Click at the bottom of the screen the  
"Queue Mode Enabled"



The "Secured System Control" Window will appear

Choose

1) Run Batch of Samples




2) Click "OK"

# Batch of Samples: Choose the "method" only one method per batch

---


The "Enter Login Name"  
Window will appear

- 1) Enter your "e-mail" ID
- 2) Enter your "Password"
- 3) Click "OK"



The "Select Analysis Method"  
Window will appear

- 1) Choose only one method
- 2) Click "OK"



*If you need to use other methods, repeat from slide #19 for each method  
(no re-login required)*

---

# Batch of Samples: Submit the Samples


---

The "MDS ChemStation" window will appear

Display the First Vial position

**Place your vial in the position indicated**

1) Click "OK"




The "Batch Sample Information" window will appear

1) Enter your "text" Sample Title

2) Enter your "text" Information

3) Enter the number of vials

4) Click "OK"



# Batch of Samples: Submit the Sample

The "MDS ChemStation" window will reiterate the number of samples – click "Yes" or "No" as appropriate

Confirm the position of the first vial, – click "OK"

Vial positions run radially from the interior to the perimeter of the vial tray




*Example:*

positions #, # + 1, # + 2

#

# + 1

# + 2



# Batch of Samples: Verify Sample Submission

After you submit the batch sample request and your sample vials have been placed in the assigned autosampler positions:

- Click on "Report Manager" and verify that your samples posted

Report Manager - "CHEMSTATION"

File Options Help

Start Pause Resume Shutdown

Posted 23 qmeth "35","Sample Name 1","Condition 1","sums" on P35667.d

# of Requests 2

System/Command	Data File	Method	Originator
RPT1 on \\CHEMSTATION\ using MSDA			
RPT2 on \\CHEMSTATION\ using MSDA2			
MC 1 on \\CHEMSTATION\ using MCTOP1			
23 qmeth "35","Sample Name 1 P35667.d		SUMS_Standard.M	CHEMSTATION
24 qmeth "36","Sample Name 2 P36668.d		SUMS_Standard.M	CHEMSTATION

Vial Positions      Filename      Method

- Important! Write the Filenames in your lab notebook



# Batch of Samples: Important reminders

---

- Never close windows on the instrument control computer! Only minimize them
  - Closing/exiting Chemstation windows will stall the queue and aggravate your colleagues
  
- All queued analyses are posted in the "Report Manager"
  - Never remove a vial in a position which is still listed in the "Report Manager"
  - If a vial position is not listed in the "Report Manager", a vial present in that position may be replaced with your sample vial


# Data Analysis

---

Open Access GC/MS

# Enhanced Data Analysis

- ❑ Windows “start” menu (if data analysis module is not already open)
- ❑ Choose “GCMS Data Analysis”
- ❑ Enter “Name” and “Password”



Enter Login Name

Name

Password

OK Cancel Help




- ❑ Click “OK”

# Open a File

- ❑ Choose File → Load Data File
- ❑ Select the File from the correct Folder  
Change Path if necessary: Choose Folder

“OK”  
and File  
“OK”





# Total Ion Chromatogram

---

- ❑ When the file is opened in Enhanced Data Analysis, the Total Ion Chromatogram (TIC) is displayed
- ❑ To zoom in on a GC peak:

Click and drag the left mouse button to draw a box around the area of interest



- ❑ A double left click on the mouse will zoom back out to the original display



# Viewing Mass Spectra

---

*Every point on the chromatogram is associated with a mass spectrum*


- To display a mass spectrum at a particular time on the TIC:
  - Position the cursor on the Chromatogram at the specific time
  - Double right click
  - The mass spectrum at that point in time will appear in a separate window
  
- To display an averaged MS across a chromatographic peak:
  - Right click and drag across the time range of interest
  
- You can toggle between windows just by clicking on them
- The active window will display a vivid blue bar at the top



# Working with Mass Spectra

---

- ❑ To zoom in: left click and drag
- ❑ To zoom out: double left click
  
- ❑ The upper mass range limit on the mass spectrum is the highest mass detected during the scan
  - Thus, even though the mass spectrometer is set to scan to the upper limit (550 Da), the x-axis of the MS often does not display up to 550.
  
- ❑ To initiate a library search:
  - Double right click on the displayed mass spectrum
  - "Wiley275" should be selected as the library; this setting is in the "select library" option in the "Spectrum" menu



# Integration

---

Open Access GC/MS




# Integration of GC peaks

---

- ❑ Chromatogram integration calculates the area under the curve of the observed GC/MS peaks
- ❑ Areas for the same compound analyzed under the same experimental condition may be related to the concentration of compound in the sample.
- ❑ *\*\*\*Appropriate controls, standards, and calibration curves should be used when quantitative measurements are desired.*
- ❑ Observed areas may be affected by sample matrix effects

# Automatic Integration

- ❑ Click on the Auto Integration button



- ❑ In the Chromatogram window, the retention time of the integrated peaks will turn blue.
- ❑ A blue line across the base of each integrated peak indicates the baseline.
- ❑ To generate a list of peak retention times and areas:
  - Select Integration Results from the Chromatogram menu
  - A Tabulated window will appear with Area of integrated peaks

# Data Reporting

---

Open Access GC/MS

# Datafile info


---

- Files are located in the folder C:\msdchem\1\DATA\“month” where “month” is the 3-letter month abbreviation.
  - Older files are in folders named “YYMM”, e.g. 0912 for Dec. 2009
- Each GC-MS data filename begins with the letter “P”; the Windows extension is “.D”
  - The data “file” is actually a folder containing the files and results for each analysis.
- Users are responsible for backing up their own data. Data files are deleted periodically as needed to free up space on the hard drive.
- Data files and data report files can be easily transferred using a flash drive (USB2)

# “Exporting” Data Views

---

- Select the window view that you would like to copy
- Click the copy button on the tool bar menu of the Enhanced Data Analysis



- Paste the clipboard on the report data file document of your choice, Office (Excel/PowerPoint/Word)


# Printing Data Views

---

- ❑ Select "Print" from the file menu
- ❑ Either or both the TIC and MS windows may be printed.
- ❑ To generate .pdf files, select "Adobe PDF" as the printer.
- ❑ When printing hard copies, ensure that the blue ethernet cable is connected (the cables to connect are located to the right of the keyboard)
- ❑ It is good practice to generate electronic data reports, rather than printing hard copies which consume paper and are easily lost.

# Exporting CSV Files

- ❑ CSV Files are comma separated value files that can be opened in Excel and saved as Excel files
- ❑ From the toolbar of the Enhanced Data Analysis window Choose: File → Export Data to CSV File
- ❑ On the Select Export Mode Window Select: Current File – With Data File
- ❑ Select the Data to export:  
Mass Spectrum  
TIC Chromatogram  
Area Percent Report (needs peak Integration)
- ❑ “OK” will generate the report data file in the datafile folder (see slide 32 for Windows file path)



# If that doesn't work...

---

- Go ask Allis.
  - [allis@stanford.edu](mailto:allis@stanford.edu)
  - 650.723.0710



Copies of this training material may be downloaded from the SUMS website at

<http://mass-spec.stanford.edu/Instruments-OpenAccess.html>

---

*Thanks to Maurizio Splendore for preparing version 1.0 of this guide, and Pavel Aronov for v 1.1.*